

Thermoplastic Splint (Nasal)

Advanced external fixation solution designed for precision, comfort, and optimal post-nasal surgery recovery.



Product Overview

Intended Use: External fixation after nose surgery.

Crafted from a special synthetic high-molecular polyester, this advanced polymer material is processed through a series of physical and chemical methods to ensure superior performance.

Its molecular structure remains stable in dry environments at room temperature (10-30°C). When heated to 65-70°C, the material softens and becomes transparent, allowing for easy stretching and precise shaping. Upon cooling, it regains rigidity while maintaining essential protective properties.

Key Characteristics



Thermoformable

Softens at 65-70°C for easy shaping, rigid at room temperature.



Breathable

Allows air circulation to maintain skin health during recovery.



Waterproof

Resistant to moisture, ensuring durability in various conditions.



Safe & Non-toxic

Odorless and non-irritating to the skin.



Biodegradable

Eco-friendly material composition.

Order Information

Shape	Size (mm)	REF Code
Trapezoid	33 × 55	Nasal-TP-01
Trapezoid	38 × 60	Nasal-TP-02
Trapezoid	45 × 60	Nasal-TP-03
Mickey Mouse	55 × 68	Nasal-MM-01
Mickey Mouse	65 × 70	Nasal-MM-02
Square	50 × 50	Nasal-SQ-01
Square	75 × 75	Nasal-SQ-02
Square	100 × 100	Nasal-SQ-03

Usage Steps



About 75°C hot water for 3 minutes.



Manufacturer & Distribution

Manufacturer

Shenzhen Alien Medical Technology Co., Ltd.

Room 1812, Building 1, Vientiane Tiancheng, Shenlong Road, Longcheng Street, Longgang District,
Shenzhen, Guangdong, China

Phone: +86 13662614220

Email: info@alienmedical.com ; tina@alienmedical.com

Website: www.alienmedical.com

EU Authorized Representative

EC REP : Riomavix Sociedad Limitada

Calle Almansa 55, 1D, Madrid 28039 Spain